

3.5 V Convertible Power Handle



71000-C Convertible Power Handle




Recyclable battery.
Please return expended rechargeable batteries to Welch Allyn for recycling.

Wherever you are,
you can power your
Welch Allyn 3.5 V convertible handle

Use either rechargeable or C-cell batteries with Welch Allyn's 3.5 V convertible power handle. Even if you've forgotten to recharge your battery or are practicing somewhere without available electricity, remove the NiCad battery and transformer module, insert two C-cell batteries and converter ring (included), and continue your diagnostic exams.

When would the convertible handle most likely be used?

The convertible handle is a lifesaver if you are on call or taking an exam and have forgotten to recharge the battery. Also, the convertible handle is very useful when you are working overseas or in areas where it is difficult to find electricity



How does this handle work with either rechargeable or C-cell batteries?

Simply remove the rechargeable battery and transformer, and replace with two C-cell batteries and the converter.



How do I recharge nickel cadmium batteries?

Our convenient built-in charger allows you to plug into any conventional 110 V wall outlet for overnight recharging.